


# Phoebus P121

## Portable Pulse Oximeter



### FEATURES

1. LCD display with in-built backlit, for easy viewing of the wave form and large display for digital read out of SpO2 and Pulse Rate value.
2. Programmable audio-visual alarm for high & low limits of SpO2 and Pulse rate.
3. Graphical trends memory of 24hrs.
4. Suitable for adult, pediatric & neonatal use.
5. Audible pitch changes with change in % of SpO2 value.
6. Mains & battery operated
7. Battery backup 4hrs.
8. Measurement Range:
  - SpO2: 0-100%
  - Pulse Rate: 30-250BPM



